

# T.E.D.<sup>™</sup> Anti-Embolism Stockings

Proven During Times of Immobility to Reduce the Likelihood of Blood Clot Formation\*



## PATIENT CONSIDERATIONS

- Are you about to have surgery?
- Are you about to go on a long flight?
- Are you likely to be immobile for long periods of time?
- Have you been identified to be at risk of a blood clot?
- Are you pregnant, on hormone replacement therapy or the contraceptive pill with one other risk factor?


If yes then the T.E.D.<sup>™</sup> Anti-Embolism Stockings may be suitable for you and reduce your risk of getting a DVT (Deep Vein Thrombosis). Ask your Doctor or Pharmacist to assess your risk.

### Indications

- To help prevent formation of emboli
- To promote increased blood flow velocity in the legs by compression of the deep venous system

### Contraindications

- Stockings are not recommended for patients with the following:
- Any local leg condition in which stockings would interfere, such as dermatitis, immediate post-operative vein ligation, gangrene or recent skin graft
  - Severe arteriosclerosis or other ischemic vascular disease
  - Massive oedema of legs or pulmonary oedema from congestive heart failure
  - Extreme deformity of leg



| KNEE LENGTH         | Calf Circumference at greatest point <b>B</b> | Heel to Bend of Knee <b>A1</b> | Knee Length Open Toe (White) | Knee Length Closed Toe (White) | Knee Length Closed Toe (Beige) | Knee Length Closed Toe (Black) |
|---------------------|---|--------------------------------|------------------------------|--------------------------------|--------------------------------|--------------------------------|
| Small Regular       | < 30.5cm                                      | 35.6cm - 40.6cm                | 7071P                        | 1046K                          | 1066                           | 1076                           |
| Small Long          | < 30.5cm                                      | > 40.6cm                       | 7339P                        | 3035                           | 3045                           | 3055                           |
| Medium Regular      | 30.5cm - 38.1cm                               | 35.6cm - 43.2cm                | 7115P                        | 1077                           | 1087K                          | 1097                           |
| Medium Long         | 30.5cm - 38.1cm                               | > 43.2cm                       | 7480P                        | 3118                           | 3128                           | 3138                           |
| Large Regular       | 38.2cm - 44.5cm                               | 35.6cm - 45.7cm                | 7203P                        | 2035                           | 2045                           | 2055                           |
| Large Long          | 38.2cm - 44.5cm                               | > 45.7cm                       | 7594P                        | 4029                           | 4039                           | 4049                           |
| Extra Large Regular | 44.6cm - 50.8cm                               | 35.6cm - 45.7cm                | 7604P                        | 2095                           | 3005                           | 3015                           |
| Extra Large Long    | 44.6cm - 50.8cm                               | > 45.7cm                       | 7802P                        | 5035                           | 5045                           | 5055                           |

| THIGH LENGTH   | Thigh Circumference at greatest point <b>C</b> | Calf Circumference at greatest point <b>B</b> | Heel to crease at buttock <b>A2</b> | Thigh Length Closed Toe |
|----------------|--|---|-------------------------------------|-------------------------|
| Small Short    | < 63.6cm                                       | < 30.5cm                                      | < 73.7cm                            | 5194                    |
| Small Regular  | < 63.6cm                                       | < 30.5cm                                      | 73.7cm - 83.8cm                     | 5213                    |
| Small Long     | < 63.6cm                                       | < 30.5cm                                      | > 83.8cm                            | 5324                    |
| Medium Short   | < 63.6cm                                       | 30.5cm - 38.1cm                               | < 73.7cm                            | 5417                    |
| Medium Regular | < 63.6cm                                       | 30.5cm - 38.1cm                               | 73.7cm - 83.8cm                     | 5533                    |
| Medium Long    | < 63.6cm                                       | 30.5cm - 38.1cm                               | > 83.8cm                            | 5628                    |
| Large Short    | < 63.6cm                                       | 38.2cm - 44.5cm                               | < 73.7cm                            | 5707                    |
| Large Regular  | < 63.6cm                                       | 38.2cm - 44.5cm                               | 73.7cm - 83.8cm                     | 5822                    |
| Large Long     | < 63.6cm                                       | 38.2cm - 44.5cm                               | > 83.8cm                            | 5943                    |

| THIGH LENGTH WITH BELT | Thigh Circumference at greatest point <b>C</b> | Calf Circumference at greatest point <b>B</b> | Heel to crease at buttock <b>A2</b> | Thigh Length (Belt) Open Toe |
|------------------------|--|---|-------------------------------------|------------------------------|
| Extra Small Regular    | < 63.6cm                                       | < 25.4cm                                      | < 71.1cm                            | 3306P                        |
| Extra Small Long       | < 63.6cm                                       | < 25.4cm                                      | 71.1cm & Longer                     | 3320P                        |
| Small Regular          | < 63.6cm                                       | 25.4cm - 30.4cm                               | < 72.5cm                            | 3039P                        |
| Small Long             | < 63.6cm                                       | 25.4cm - 30.4cm                               | 72.5cm & Longer                     | 3364P                        |
| Medium Regular         | < 63.6cm                                       | 30.5cm - 38.1cm                               | < 72.5cm                            | 3144P                        |
| Medium Long            | < 63.6cm                                       | 30.5cm - 38.1cm                               | 72.5cm & Longer                     | 3449P                        |
| Large Regular          | < 63.6cm                                       | 38.2cm - 44.5cm                               | < 73.8cm                            | 3221P                        |
| Large Long             | < 63.6cm                                       | 38.2cm - 44.5cm                               | 73.8cm & Longer                     | 3523P                        |
| Extra Large Regular    | 63.6cm - 81.3cm                                | 38.2cm - 44.5cm                               | < 72.5cm                            | 3922P                        |
| Extra Large Long       | 63.6cm - 81.3cm                                | 38.2cm - 44.5cm                               | 72.5cm & Longer                     | 3995P                        |


## HOW TO MEASURE CORRECTLY

### KNEE LENGTH MEASUREMENTS

1. Measure calf circumference at greatest portion to determine size. **B**
2. Measure the distance from bend of knee to bottom of heel to determine length. If possible, measure length in standing position. **A1**

### THIGH LENGTH MEASUREMENTS

1. Measure upper thigh circumference at the buttock fold. When thigh circumference measures greater than 91.4 cm, select knee length style. **C**
2. Measure calf circumference at greatest portion to determine size. **B**
3. Measure distance from base of heel to gluteal furrow. **A2**



1: Wells, PS; Lensing, AWA; Hirsh J. Graduated Compression Stockings in the Prevention of Postoperative Venous Thromboembolism. A Meta-analysis. Arch. Intern. Med. 1994; 154: 67-71.